

WIRELESS IN-BORE VIEWING



Head-First Scan - MRI Wireless Portable Video Display (WPVD) is placed behind the bore. Patient uses mirror glasses, head-coil mirror or bridge mirror to see display

“They aren’t thinking about going in the tube, they are thinking about...the color, video and music. They like all of those different options.”

- Rachel, MRI Tech



Feet-First Scan - MRI Wireless Portable Video Display can be placed next to the patient or behind the bore

MRI-SAFE WIRELESS BATTERY POWERED PORTABLE VIDEO DISPLAY (MRI WPVD)

- MRI-safe, no cords or cables
- Easily moves around the scan room
- Operates for 13 hours on full battery charge
- Recharge MRI WPVD overnight or operate with charger plugged in
- Compatible with all MRI systems
- Audio connects to the MRI patient audio system



PDC MRI WPVD



MRI WPVD PROVEN BENEFITS:

- Reduces sedation and anesthesia
- Improves MRI through-put
- Increases MRI utilization and revenue
- Improves technologist workflow
- Improves patient satisfaction
- Integrates with fMRI

TWO MRI WPVD SOLUTIONS AVAILABLE:

Essential MRI WPVD Solution -

- Easy three-step installation
- Connects to customer's audio video source, patient audio system, smartphone, blu-ray or any streaming device
- Installation complete in less than 2 hours
- No construction required

Premium MRI WPVD Solution -

- Includes PDC's control system
- Touchpad control
- Streaming capability
- Original PDC Nature Theme videos
- Custom scrolling text on MRI WPVD
- Countdown timer on MRI WPVD
- Installation requires mounting control system in the equipment room
- Optional scan room speakers




THE ESSENTIAL INSTALLATION IS COMPLETE IN THREE EASY STEPS - NO CONSTRUCTION REQUIRED

Control Room Glass Wireless Transmitter (Tx)



1. Place Tx in front of RF glass in control room. Plug Tx into a standard wall outlet

Scan Room Wireless Transmitter/Receiver (TxRx)



2. Place TxRx in the scan room against RF glass opposite the Tx, plug into a standard wall outlet

MRI Wireless PVD



3. Position MRI WPVD in scan room behind magnet and press power button