

TO POWER WIRELESS IN-BORE VIEWING VIDEO DISPLAY (WIRELESS VIDEO DISPLAY) ON/OFF:

- Press power button on top of Wireless Video Display to power ON and OFF
- Power button illuminates when powered ON
- Power OFF and charging recommended nightly

TO CHARGE:


1. Plug charging cable into top of Wireless Video Display until it clicks
2. Plug other end of charging cable into Wireless Video Display outlet
3. Charging can take place during an MRI scan
4. With the Wireless Video Display powered OFF, a full charge takes about 8 hours
5. Charging nightly with Wireless Video Display powered OFF is recommended
6. See indicators below for charge status

WIRELESS LED INDICATORS:

- 3 Lit **GREEN** LEDs = Fully charged & charger is connected
- 2 Lit **GREEN** LEDs = Fully charged & charger is disconnected. ~12 hrs battery life
- 1 Lit **GREEN** LED = Half charge. Approximately 6 hours of battery life remaining.
- 1 Lit **RED** LED = Battery is depleted and requires recharge. Monitor will power OFF when battery reaches this point

SCAN FOR WIRELESS VIDEO DISPLAY
TRAINING & MARKETING RESOURCES





POSITIONING:

- When possible, keep Wireless Video Display 3 feet from magnet
- For in-bore viewing, position Wireless Video Display behind magnet bore, at least 3 feet from bore
- To adjust the height of Wireless Video Display, turn the **BLACK** lever on the stand, closest to display, to the noted unlocked position, adjust height and return lever to locked position
- Lock wheels in place by pressing "down" on wheel levers. Lift levers "up" to unlock

BATTERY CHARGER:

- For MRI, located in equipment room. For other modalities, can be placed anywhere
- When Wireless Video Display is connected to charger, charger will light up **RED** to indicate charging
- Ensure charger is set to **10A** charge rate

TROUBLESHOOTING:

- If Wireless Video Display says "Waiting to Connect," move Wireless Video Display away from magnet and power OFF and ON using the Wireless Video Display power button
- If Wireless Video Display says, "No input," check to ensure a source is selected on the Caring Suite™ App and consult App troubleshooting if needed

SCAN FOR MRI CARING SUITE
TRAINING & MARKETING RESOURCES

