PDC

LIGHTING + AUDIO CARING SUITE™

=

One-Touch Controls		Suite Levels	0 9 0
0		Main Color	Power Options
0	E	Ceiling Image	Refresh App
		Can Lights	Reload App O Restart Server
<		-	(b) Restart Tablet Toolkit
		 Directional Lights - IIII 	(c) Reboot Tablet
Start Scan ighting Suite	Patient Prep		Power Down Suite

TROUBLESHOOTING

- FOR CONNECTION ISSUES/RED RACK ICON: 1st, tap Power options at bottom right & tap "Restart Tablet Toolkit" or "Reboot Tablet"
- FOR GENERAL APP ISSUES: 1st tap "Refresh App." If not successful, tap "Reload App." If issue is not resolved, tap "Restart Server"
- FOR ADDITIONAL SUPPORT contact PDC Tech Support, info below

