



## IN-BORE PROJECTOR - SIEMENS WALL MOUNT

- **POSITIONING THE MIRROR LENS ASSEMBLY (MLA):** Secure on table as with Siemens' mirror. For use with head coil, slide fully back to a hard stop. If patient's head is not positioned near table end, slide forward until image is in front of patient's face. Remove MLA before lowering table
- **POWERING OFF / ON PROJECTOR:** To power on / off Projector, press the power button on the bottom side of the projector in image 1, or unplug, if needed. Power off recommended nightly, when not in use - or once per day
- **SELECTING PROJECTOR VIDEO:** On the Caring Suite™ App, tap the blue "Video Control" box. Next, tap the blue bar below "In-Bore Display" (name could vary). Tap to select the desired source, i.e., "Themes" or "DVD w/Sound." See p11 for more
- **TROUBLESHOOTING:** If the Projector image says "Waiting for Connection" or is not displaying video: 1) Power off Projector. 2) Power off system via App Power button at lower right. Power back on system and power back on Projector. If not successful, power off Projector and follow troubleshooting steps for Power Events. Once complete, power on Projector.



## TROUBLESHOOTING: POWER EVENTS

### • FOR PLANNED OUTAGES:

- If possible, locate System Control Rack in equipment area
- Power OFF the computer by pressing the white power button to right of the red PDC. Ensure button is no longer illuminated, signaling the computer is OFF
- Flip power switch to left of red PDC to power OFF audio / video
- Once power surge or outage has passed, flip power switch back to ON position

### • IF AUDIO / VIDEO DOES NOT RETURN AFTER A POWER EVENT

- In App power options at bottom right, tap "Restart Server." If not successful:
- Locate \*System Control Rack in equipment area
- Power OFF computer by pressing the white power button to right of red PDC
- Flip power switch to left of PDC to OFF. Wait ~30 seconds & power back on
- Power cycle AVB (Audio Video Box) below tablet by unplugging power cable at the back, waiting about a minute and replugging the power cable

### • IF COMPUTER WON'T POWER ON:

- Check the CyberPower UPS on bottom or side of rack. A blue light will indicate it is ON. If it is OFF, open the window covering the button & press to turn it ON. If it's ON, check to ensure it's sending power by pressing the silver button below the power button. Check to ensure "Input" & "Output" power is ~120v ac

\*Note: System Control Rack may be divided into 2 parts