## $\bigcirc \mathbf{PDC} \quad \mathsf{CARING SUITE^{TM} APP GUIDE}$

	নি: [	🖗 One-Touch C	onfiguration				$\rightarrow$	?	۲	
		C Restore Configuration Build Default Remove a Control Save Configuration								
		One-Touch Controls			One-Touch Options					
		<del>مۇ</del>	<ul><li>✓ ~€</li></ul>	A/V Routes	GON-Screen Dis	play Profiles				
s	Suit	ViewRay M2 Standal (Portable Displays)	M2 Standard ViewRay M2 Zoom le Displays) (Portable Displays)	o to					ł	
		*	-t:	✓	Themes (Suite Audio)					
		ViewRay M2 RealVie (Portable Displays)	w Themes ) (Suite Audio)	-to						
		•		P (P	ortable Camera ortable Displays)					
Admin Portal		5 Seconds	Please hold still	1	$\bigcirc$					

## CUSTOMIZING ONE-TOUCH CONTROLS

- TAP MENU icon in upper left & tap "Administrator Portal." Enter: "8675309" when prompted to login
- TAP "APP CONFIGURATION." Below Tech Portal, tap "Configure" beside "One-Touch Configuration"
- TAP "REMOVE A CONTROL" at top, & tap to remove an unwanted option from the left below "One-Touch Controls." Repeat for all unwanted options
- FIND DESIRED OPTION(S) at right. Available options are noted with a "+" symbol to the right. Tap the "+" symbol beside the desired option to add
- TO SAVE CHANGES tap "Save Configuration" at top, tap Menu icon at left & tap "Tech Portal" to return to main view
- TO RESTORE DEFAULT tap "Build Default" at top

