



MIRROR OPTIONS FOR INBORE VIEWING

- **MIRROR PLATFORM:** Designed to keep mirror stable on patient table
- **PROTECTIVE PLASTIC:** To be removed before using the mirror
- **ADJUSTABLE ARM:** Loosen black knob to adjust, if needed
- **ADJUSTABLE MIRROR:** The patient can adjust the mirror to optimize their view after they are in position on the table
- **MIRROR POSITIONING:** Place the platform on the patient table where the patient's head will rest & cover with a sheet, pillow or soft pad so the patient is comfortable
- **PLATFORM PEGS:** Fit into indentations on the edge of the MRIdian patient table to help keep the mirror stable
- **SPECIFIC TO MRIDIAN®:** This mirror was developed specifically for the MRIdian
- **MIRROR GLASSES:** Mirror can be manually tilted to adjust view