



PDC CARING SUITETM APP GUIDE




ViewRay M2 - Mirrors ViewRay Computer



ViewRay M2 Standard - For Treatment



ViewRay M2 Zoom - For Treatment



ViewRay RealView - For Simulation

VIDEO CONTROL + TREATMENT VIEWS

- YOUR CARING SUITE INCLUDES 4 VIEWRAY OPTIONS:
 - VIEWRAY M2 reflects your ViewRay computer screen
 - VIEWRAY M2 STANDARD shows a custom view for treatment, as above
 - VIEWRAY M2 ZOOM shows a custom zoomed view for treatment, as above
 - VIEWRAY REALVIEW shows a custom view for use during simulation, as above
- TO SELECT VIEWRAY OPTIONS: See P10
- TO ADD OR REMOVE OPTIONS FROM ONE-TOUCH CONTROLS: See P6

